

STEREOMICROSCOPES with CAMERA

The following illustrations show the model SMZ zoom 143 with the triokulartubus. The camera Moticam 1000 was attached. Fotos and films can be transferred in high resolution to computers and screens. Optimal for projector 1280 x 1024 pixels. The provided software for Windows and Apple Mac OS X makes beyond that over the standardized USB connection versatile measurements possible.



Stereomicroscop SMZ Zoom 143 with Camera Moticam 1000

The connection between camera and optics is made with CCD adapter 0,5x, which is provided as accessories to the microscope.

Art.-No. 2.10


The camera **Moticam 1000** is delivered with several eyepiece adapters in the dimension 28, 30 and 34mm so that the camera can be fastened also without okulartubus directly to the eyepiece on usual standard microscopes.

Art. - No. 2.20



Stereomicroscop SMZ Zoom 143 with Moticam 1000

Attached is the Apple-computer iMac. On the screen the pressure grip forceps is shown. Under the optics is the instrument 1.02.



Successor is the new model **Motic 1 SP**.

The presented articles are specified in the price list.

